

MLKL Rabbit mAb

Catalog No: #60050

Package Size: #60050-1 50ul #60050-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MLKL Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	MLKL Antibody detects endogenous levels of total MLKL
Immunogen Description	A synthesized peptide derived from MLKL
Other Names	hMLKL; Mixed lineage kinase domain like; Mlk1;
Accession No.	Uniprot:Q8NB16
Uniprot	Q8NB16
Calculated MW	54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details


WB 1:500~1:2000 IP 1:50

Images



Western blot analysis of MLKL expression in HeLa cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Product Description

Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process.

Background

Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process.

Note: This product is for in vitro research use only