

UBE4B Rabbit mAb

Catalog No: #60061

Package Size: #60061-1 50ul #60061-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	UBE4B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	UBE4B Antibody detects endogenous levels of total UBE4B
Immunogen Description	A synthesized peptide derived from UBE4B
Other Names	HDNB1; Ube4b; UBOX3; UFD2; UFD2a; Ufd2p;
Accession No.	Uniprot:O95155
Uniprot	O95155
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of UBE4B expression in (1) 293T cell lysate; (2) NIH/3T3 cell lysate.

Product Description

Binds to the ubiquitin moieties of preformed conjugates and catalyzes ubiquitin chain assembly in conjunction with E1, E2, and E3.

Background

Binds to the ubiquitin moieties of preformed conjugates and catalyzes ubiquitin chain assembly in conjunction with E1, E2, and E3.

Note: This product is for in vitro research use only