

## PGD Rabbit mAb

Catalog No: #60065

Package Size: #60065-1 50ul #60065-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PGD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	PGD Antibody detects endogenous levels of total PGD
Immunogen Description	A synthesized peptide derived from PGD
Other Names	6PGD; Pgd;
Accession No.	Uniprot:P52209
Uniprot	P52209
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

## Images



Western blot analysis of PGD expression in (1) 293 cell lysate;  
(2) Mouse kidney lysate; (3) Rat kidney lysate.

## Product Description

Catalyzes the oxidative decarboxylation of 6-phosphogluconate to ribulose 5-phosphate and CO<sub>2</sub>, with concomitant reduction of NADP to NADPH.

## Background

Catalyzes the oxidative decarboxylation of 6-phosphogluconate to ribulose 5-phosphate and CO<sub>2</sub>, with concomitant reduction of NADP to NADPH.

---

Note: This product is for in vitro research use only