

MGST1 Rabbit mAb

Catalog No: #60071

Package Size: #60071-1 50ul #60071-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MGST1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	MGST1 Antibody detects endogenous levels of total MGST1
Immunogen Description	A synthesized peptide derived from human MGST1
Other Names	MGST1; GST12; MGST; MGST I; MGST1;
Accession No.	Uniprot:P10620
Uniprot	P10620
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of MGST1 expression in HepG2 cell lysate.

Product Description

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Has a wide substrate specificity.

Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Has a wide substrate specificity.

Note: This product is for in vitro research use only