

## C1QBP Rabbit mAb

Catalog No: #60073

Package Size: #60073-1 50ul #60073-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	C1QBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	C1QBP Antibody detects endogenous levels of total C1QBP
Immunogen Description	A synthesized peptide derived from human C1QBP
Other Names	C1qBP; GC1QBP; GC1QR; HABP1; p32 ; p33; SF2p32;
Accession No.	Uniprot:Q07021
Uniprot	Null
Calculated MW	31kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details


WB 1:500~1:2000 IHC 1:50~1:200

## Images



Western blot analysis of C1QBP expression in HeLa cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



## Product Description

Binds to the globular heads" of C1Q thus inhibiting C1 activation.....1,,1,1,1,0,1,1701055481,0,0,0,0,0,0,1701055481,0,-1,-1,0,0  
324741,57,60074,ORP1 Rabbit mAb,0,10,ORP1-Rabbit-mAb,ORP1 Rabbit  
mAb,,0,50ul,249,249,0,0,0,58,ABWAY,CY10417,CY10417,0,0,0,Yes,Rabbit Monoclonal Antibody,Binds phospholipids; exhibits strong binding to  
phosphatidic acid and weak binding to phosphatidylinositol 3-phosphate.....1,,1,1,1,0,1,1701055481,0,0,0,0,0,0,1701055481,0,-1,-1,0,0  
324742,57,60075,MUTA Rabbit mAb,0,10,MUTA-Rabbit-mAb,MUTA Rabbit  
mAb,,0,50ul,249,249,0,0,0,58,ABWAY,CY10441,CY10441,0,0,0,Yes,Rabbit Monoclonal Antibody,"Involved in the degradation of several amino acids,  
odd-chain fatty acids and cholesterol via propionyl-CoA to the tricarboxylic acid cycle. MCM has different functions in other species.

## Background

Binds to the globular heads" of C1Q thus inhibiting C1 activation.,,Q07021,Null,<http://www.uniprot.org/uniprot/Q07021>,Null  
324741,324741,60074,Null,ORP1 Rabbit mAb,Rabbit,Monoclonal,Null,Rabbit IgG,Affinity-chromatography,WB ICC/IF FC,WB 1:500~1:2000 ICC/IF  
1:50~1:200 FC 1:50,Human Mouse Rat ,ORP1 Antibody detects endogenous levels of total ORP1,,,,,ORP1; OSBP8; OSBPL1; OSBPL1A;  
OSBPL1B;Uniprot:Q9BXW6,,Null,Null,"50,108kDa

## References

Null

Note: This product is for in vitro research use only