

SCIN Rabbit mAb

Catalog No: #60081

Package Size: #60081-1 50ul #60081-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SCIN Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	SCIN Antibody detects endogenous levels of total SCIN
Immunogen Description	A synthesized peptide derived from human SCIN
Other Names	Adseverin; Scin; Scinderin;
Accession No.	Uniprot:Q9Y6U3
Uniprot	Q9Y6U3
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of SCIN expression in mouse kidney cell lysate.

Product Description

Ca²⁺-dependent actin filament-severing protein that is presumed to have a regulatory function in exocytosis by affecting the organization of the microfilament network underneath the plasma membrane.

Background

Ca(2+)-dependent actin filament-severing protein that is presumed to have a regulatory function in exocytosis by affecting the organization of the microfilament network underneath the plasma membrane.

Note: This product is for in vitro research use only