

## Lactate Dehydrogenase C Rabbit mAb

Catalog No: #60089

Package Size: #60089-1 50ul #60089-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Lactate Dehydrogenase C Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Lactate Dehydrogenase C Antibody detects endogenous levels of total Lactate Dehydrogenase C
Immunogen Description	A synthesized peptide derived from human Lactate Dehydrogenase C
Other Names	CT32; Lactate dehydrogenase C4; LDH3; Idhc; LDHX;
Accession No.	Uniprot:P07864
Uniprot	P07864
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of Lactate Dehydrogenase C expression in human testis lysate.

## Product Description

Possible role in sperm motility.

## Background

Possible role in sperm motility.

---

Note: This product is for in vitro research use only