

HMG4 Rabbit mAb

Catalog No: #60114

Package Size: #60114-1 50ul #60114-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|--|
| Product Name | HMG4 Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Purification | Affinity-chromatography |
| Applications | WB |
| Species Reactivity | Human Mouse |
| Specificity | HMG4 Antibody detects endogenous levels of total HMG4 |
| Immunogen Description | A synthesized peptide derived from human HMG4 |
| Other Names | HMG 4; HMG2A; Hmgb3; |
| Accession No. | Uniprot:O15347 |
| Uniprot | O15347 |
| Formulation | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle. |

Application Details

WB 1:500~1:2000

Images



Western blot analysis of HMG4 expression in HeLa cell lysate.

Product Description

Binds preferentially single-stranded DNA and unwinds double stranded DNA.

Background

Binds preferentially single-stranded DNA and unwinds double stranded DNA.

Note: This product is for in vitro research use only