

## CYB5R3 Rabbit mAb

Catalog No: #60118

Package Size: #60118-1 50ul #60118-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CYB5R3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Rat
Specificity	CYB5R3 Antibody detects endogenous levels of total CYB5R3
Immunogen Description	A synthesized peptide derived from human CYB5R3
Other Names	B5R; Cyb5r3; Cytochrome b5 reductase 3; DIA1; Diaphorase 1;
Accession No.	Uniprot:P00387
Uniprot	P00387
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200

## Images



Western blot analysis of CYB5R3 expression in A431 cell lysate.

## Product Description

Desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction.

## Background

Desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction.

---

Note: This product is for in vitro research use only