

CYB5R3 Rabbit mAb

Catalog No: #60118

Package Size: #60118-1 50ul #60118-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CYB5R3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Rat
Specificity	CYB5R3 Antibody detects endogenous levels of total CYB5R3
Immunogen Description	A synthesized peptide derived from human CYB5R3
Other Names	B5R; Cyb5r3; Cytochrome b5 reductase 3; DIA1; Diaphorase 1;
Accession No.	Uniprot:P00387
Uniprot	P00387
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of CYB5R3 expression in A431 cell lysate.

Product Description

Desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction.

Background

Desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction.

Note: This product is for in vitro research use only