

TORC2 Rabbit mAb

Catalog No: #60119

Package Size: #60119-1 50ul #60119-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TORC2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	TORC2 Antibody detects endogenous levels of total TORC2
Immunogen Description	A synthesized peptide derived from human TORC2
Other Names	CRTC2; torc2; Transducer of CREB protein 2;
Accession No.	Uniprot:Q53ET0
Uniprot	Q53ET0
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of TORC2 expression in 293T cell lysate.

Product Description

Transcriptional coactivator for CREB1 which activates transcription through both consensus and variant cAMP response element (CRE) sites.

Background

Transcriptional coactivator for CREB1 which activates transcription through both consensus and variant cAMP response element (CRE) sites.

Note: This product is for in vitro research use only