

HDLBP Rabbit mAb

Catalog No: #60121

Package Size: #60121-1 50ul #60121-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HDLBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	HDLBP Antibody detects endogenous levels of total HDLBP
Immunogen Description	A synthesized peptide derived from human HDLBP
Other Names	HBP; HDL binding protein; HDLBP; PRO2900; VGL; Vigilin;
Accession No.	Uniprot:Q00341
Uniprot	Q00341
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of HDLBP expression in HeLa cell lysate.

Product Description

Appears to play a role in cell sterol metabolism. It may function to protect cells from over-accumulation of cholesterol.

Background

Appears to play a role in cell sterol metabolism. It may function to protect cells from over-accumulation of cholesterol.

Note: This product is for in vitro research use only