

Cadherin 10 Rabbit mAb

Catalog No: #60123

Package Size: #60123-1 50ul #60123-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Cadherin 10 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Cadherin 10 Antibody detects endogenous levels of total Cadherin 10
Immunogen Description	A synthesized peptide derived from human Cadherin 10
Other Names	Cadherin 10 type 2; CDH10; T2 cadherin;
Accession No.	Uniprot:Q9Y6N8
Uniprot	Q9Y6N8
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Cadherin 10 expression in U-87 MG cell lysate.

Product Description

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Background

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Note: This product is for in vitro research use only