

## Pleckstrin Rabbit mAb

Catalog No: #60130

Package Size: #60130-1 50ul #60130-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Pleckstrin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Pleckstrin Antibody detects endogenous levels of total Pleckstrin
Immunogen Description	A synthesized peptide derived from human Pleckstrin
Other Names	p47; Pleckstrin; PLEK;
Accession No.	Uniprot:P08567
Uniprot	P08567
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of Pleckstrin expression in Raji cell lysate.

## Product Description

Major protein kinase C substrate of platelets.

## Background

Major protein kinase C substrate of platelets.

---

Note: This product is for in vitro research use only