

BCKDK Rabbit mAb

Catalog No: #60133

Package Size: #60133-1 50ul #60133-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BCKDK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	BCKDK Antibody detects endogenous levels of total BCKDK
Immunogen Description	A synthesized peptide derived from human BCKDK
Other Names	BCKD kinase; BCKDHKIN; Bckdk; BCKDKD; BDK;
Accession No.	Uniprot:O14874
Uniprot	O14874
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of BCKDK expression in HeLa cell lysate.

Product Description

Catalyzes the phosphorylation and inactivation of the branched-chain alpha-ketoacid dehydrogenase complex, the key regulatory enzyme of the valine, leucine and isoleucine catabolic pathways. Key enzyme that regulate the activity state of the BCKD complex.

Background

Catalyzes the phosphorylation and inactivation of the branched-chain alpha-ketoacid dehydrogenase complex, the key regulatory enzyme of the valine, leucine and isoleucine catabolic pathways. Key enzyme that regulate the activity state of the BCKD complex.

Note: This product is for in vitro research use only