

SSTR3 Rabbit mAb

Catalog No: #60140

Package Size: #60140-1 50ul #60140-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SSTR3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human
Specificity	SSTR3 Antibody detects endogenous levels of total SSTR3
Immunogen Description	A synthesized peptide derived from human SSTR3
Other Names	Smstr3; SS3R; SSR28; Sst3; SSTR3;
Accession No.	Uniprot:P32745
Uniprot	P32745
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 FC 1:50

Images



Western blot analysis of SSTR3 expression in Somatostatin Receptor 3 SSTR3 transfected 293T cell lysate.

Product Description

Receptor for somatostatins-14 and -28. This receptor is coupled via pertussis toxin sensitive G proteins to inhibition of adenylyl cyclase.

Background

Receptor for somatostatins-14 and -28. This receptor is coupled via pertussis toxin sensitive G proteins to inhibition of adenylyl cyclase.

Note: This product is for in vitro research use only