

MRPS31 Rabbit mAb

Catalog No: #60153

Package Size: #60153-1 50ul #60153-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MRPS31 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	MRPS31 Antibody detects endogenous levels of total MRPS31
Immunogen Description	A synthesized peptide derived from human MRPS31
Other Names	IMOGN38; MRPS31; S31mt;
Accession No.	Uniprot:Q92665
Uniprot	Q92665
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of MRPS31 expression in HepG2 cell lysate.

Product Description

Component of the mitochondrial ribosome small subunit (28S) which comprises a 12S rRNA and about 30 distinct proteins.

Background

Component of the mitochondrial ribosome small subunit (28S) which comprises a 12S rRNA and about 30 distinct proteins.

Note: This product is for in vitro research use only