

WDR1 Rabbit mAb

Catalog No: #60159

Package Size: #60159-1 50ul #60159-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	WDR1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	WDR1 Antibody detects endogenous levels of total WDR1
Immunogen Description	A synthesized peptide derived from human WDR1
Other Names	AIP1; NORI 1; wdr1;
Accession No.	Uniprot:O75083
Uniprot	O75083
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

Images



Western blot analysis of WDR1 expression in Raji cell lysate.

Product Description

Induces disassembly of actin filaments in conjunction with ADF/cofilin family proteins (PubMed:15629458). Enhances cofilin-mediated actin severing (By similarity).

Background

Induces disassembly of actin filaments in conjunction with ADF/cofilin family proteins (PubMed:15629458). Enhances cofilin-mediated actin severing (By similarity).

Note: This product is for in vitro research use only