

## NgR3 Rabbit mAb

Catalog No: #60165

Package Size: #60165-1 50ul #60165-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	NgR3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	NgR3 Antibody detects endogenous levels of total NgR3
Immunogen Description	A synthesized peptide derived from human NgR3
Other Names	NgR3; NGRH2; NGRL2; RTN4RL1;
Accession No.	Uniprot:Q86UN2
Uniprot	Q86UN2
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of NgR3 expression in K562 cell lysate.

## Product Description

May play a role in regulating axonal regeneration and plasticity in the adult central nervous system.

## Background

May play a role in regulating axonal regeneration and plasticity in the adult central nervous system.

---

Note: This product is for in vitro research use only