

DLX3 Rabbit mAb

Catalog No: #60166

Package Size: #60166-1 50ul #60166-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DLX3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	DLX3 Antibody detects endogenous levels of total DLX3
Immunogen Description	A synthesized peptide derived from human DLX3
Other Names	AI4; Dlx3; TDO;
Accession No.	Uniprot:O60479
Uniprot	O60479
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of DLX3 expression in JAR cell lysate.

Product Description

Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis.

Background

Likely to play a regulatory role in the development of the ventral forebrain. May play a role in craniofacial patterning and morphogenesis.

Note: This product is for in vitro research use only