

PEX11B Rabbit mAb

Catalog No: #60168

Package Size: #60168-1 50ul #60168-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PEX11B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	PEX11B Antibody detects endogenous levels of total PEX11B
Immunogen Description	A synthesized peptide derived from human PEX11B
Other Names	Peroxin 11B; PEX11 beta; PEX11B; Pex11p beta; PEX14B;
Accession No.	Uniprot:O96011
Uniprot	O96011
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of PEX11B expression in HT-29 cell lysate.

Product Description

Involved in peroxisomal proliferation. May regulate peroxisomes division by recruiting the dynamin-related GTPase DNM1L to the peroxisomal membrane.

Background

Involved in peroxisomal proliferation. May regulate peroxisomes division by recruiting the dynamin-related GTPase DNM1L to the peroxisomal membrane.

Note: This product is for in vitro research use only