

Nup107 Rabbit mAb

Catalog No: #60169

Package Size: #60169-1 50ul #60169-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Nup107 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Nup107 Antibody detects endogenous levels of total Nup107
Immunogen Description	A synthesized peptide derived from human Nup107
Other Names	NUP107; NUP84;
Accession No.	Uniprot:P57740
Uniprot	P57740
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of Nup107 expression in 293 cell lysate.

Product Description

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Required for the assembly of peripheral proteins into the NPC. May anchor NUP62 to the NPC.

Background

Plays a role in the nuclear pore complex (NPC) assembly and/or maintenance. Required for the assembly of peripheral proteins into the NPC. May anchor NUP62 to the NPC.

Note: This product is for in vitro research use only