

TCEA1 Rabbit mAb

Catalog No: #60182

Package Size: #60182-1 50ul #60182-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCEA1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	TCEA1 Antibody detects endogenous levels of total TCEA1
Immunogen Description	A synthesized peptide derived from human TCEA1
Other Names	GTF2S; TCEA; Tcea1; TF2S; TFIIS;
Accession No.	Uniprot:P23193
Uniprot	P23193
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of TCEA1 expression in HeLa cell lysate.

Product Description

Necessary for efficient RNA polymerase II transcription elongation past template-encoded arresting sites.

Background

Necessary for efficient RNA polymerase II transcription elongation past template-encoded arresting sites.

Note: This product is for in vitro research use only