

## TEMT Rabbit mAb

Catalog No: #60183

Package Size: #60183-1 50ul #60183-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TEMT Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	TEMT Antibody detects endogenous levels of total TEMT
Immunogen Description	A synthesized peptide derived from human TEMT
Other Names	Amine N methyltransferase; Arylamine N methyltransferase; Indolethylamine N methyltransferase; Inmt; TEMT; Thioether S methyltransferase;
Accession No.	Uniprot:O95050
Uniprot	O95050
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of TEMT expression in A549 cell lysate.

## Product Description

Catalyzes the N-methylation of tryptamine and structurally related compounds.

## Background

---

Catalyzes the N-methylation of tryptamine and structurally related compounds.

---

Note: This product is for in vitro research use only