

FPGS Rabbit mAb

Catalog No: #60184

Package Size: #60184-1 50ul #60184-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FPGS Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	FPGS Antibody detects endogenous levels of total FPGS
Immunogen Description	A synthesized peptide derived from human FPGS
Other Names	FOLC; Folylpoly gamma glutamate synthetase; FPGS; Tetrahydrofolate synthase;
Accession No.	Uniprot:Q05932
Uniprot	Q05932
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of FPGS expression in HeLa cell lysate.

Product Description

Conversion of folates to polyglutamate derivatives. This allows tissues to concentrate folate at higher levels than in plasma.

Background

Conversion of folates to polyglutamate derivatives. This allows tissues to concentrate folate at higher levels than in plasma.

Note: This product is for in vitro research use only