

## IFT88 Rabbit mAb

Catalog No: #60185

Package Size: #60185-1 50ul #60185-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	IFT88 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	IFT88 Antibody detects endogenous levels of total IFT88
Immunogen Description	A synthesized peptide derived from human IFT88
Other Names	DAF19; hTg737; Ift88; TG737; TTC10;
Accession No.	Uniprot:Q13099
Uniprot	Q13099
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of IFT88 expression in HepG2 cell lysate.

## Product Description

Involved in primary cilium biogenesis.

## Background

Involved in primary cilium biogenesis.

---

Note: This product is for in vitro research use only