

## PRPF31 Rabbit mAb

Catalog No: #60189

Package Size: #60189-1 50ul #60189-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PRPF31 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	PRPF31 Antibody detects endogenous levels of total PRPF31
Immunogen Description	A synthesized peptide derived from human PRPF31
Other Names	hPrp31; Protein 61K; PRP31; prpf31; RP11; SNRNP61;
Accession No.	Uniprot:Q8WWY3
Uniprot	Q8WWY3
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

## Images



Western blot analysis of PRPF31 expression in Jurkat cell lysate.

## Product Description

Involved in pre-mRNA splicing. Required for U4/U6.U5 tri-snRNP formation.

## Background

Involved in pre-mRNA splicing. Required for U4/U6.U5 tri-snRNP formation.

---

Note: This product is for in vitro research use only