

MSRB3 Rabbit mAb

Catalog No: #60190

Package Size: #60190-1 50ul #60190-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MSRB3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	MSRB3 Antibody detects endogenous levels of total MSRB3
Immunogen Description	A synthesized peptide derived from human MSRB3
Other Names	DFNB74; MsrB3;
Accession No.	Uniprot:Q8IXL7
Uniprot	Q8IXL7
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of MSRB3 expression in HeLa cell lysate.

Product Description

Catalyzes the reduction of free and protein-bound methionine sulfoxide to methionine. Isoform 2 is essential for hearing.

Background

Catalyzes the reduction of free and protein-bound methionine sulfoxide to methionine. Isoform 2 is essential for hearing.

Note: This product is for in vitro research use only