

TMEM192 Rabbit mAb

Catalog No: #60194

Package Size: #60194-1 50ul #60194-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TMEM192 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	TMEM192 Antibody detects endogenous levels of total TMEM192
Immunogen Description	A synthesized peptide derived from human TMEM192
Other Names	TMEM192; Transmembrane protein 192;
Accession No.	Uniprot:Q8IY95
Uniprot	Q8IY95
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of TMEM192 expression in U87-MG cell lysate.

Product Description

Belongs to the TMEM192 family.

Background

Belongs to the TMEM192 family.

Note: This product is for in vitro research use only