

MYO1B Rabbit mAb

Catalog No: #60201

Package Size: #60201-1 50ul #60201-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MYO1B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	MYO1B Antibody detects endogenous levels of total MYO1B
Immunogen Description	A synthesized peptide derived from human MYO1B
Other Names	MMI alpha; MMIa; MYH 1c; Myo1b; myr1;
Accession No.	Uniprot:O43795
Uniprot	O43795
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of MYO1B expression in (1) 293 cell lysate; (2) Mouse heart lysate.

Product Description

Motor protein that may participate in process critical to neuronal development and function such as cell migration, neurite outgrowth and vesicular transport.

Background

Motor protein that may participate in process critical to neuronal development and function such as cell migration, neurite outgrowth and vesicular transport.

Note: This product is for in vitro research use only