

TNS4 Rabbit mAb

Catalog No: #60204

Package Size: #60204-1 50ul #60204-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TNS4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	TNS4 Antibody detects endogenous levels of total TNS4
Immunogen Description	A synthesized peptide derived from human TNS4
Other Names	CTEN; Tensin 4; Tns4;
Accession No.	Uniprot:Q8IZW8
Uniprot	Q8IZW8
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of TNS4 expression in HepG2 cell lysate.

Product Description

May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton (By similarity). May promote apoptosis, via its cleavage by caspase-3.

Background

May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton (By similarity). May promote apoptosis, via its cleavage by caspase-3.

Note: This product is for in vitro research use only