

## RNF115 Rabbit mAb

Catalog No: #60208

Package Size: #60208-1 50ul #60208-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	RNF115 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	RNF115 Antibody detects endogenous levels of total RNF115
Immunogen Description	A synthesized peptide derived from human RNF115
Other Names	BCA2; RNF115;
Accession No.	Uniprot:Q9Y4L5
Uniprot	Q9Y4L5
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of RNF115 expression in PC-3 cell lysate.

## Product Description

Acts as an E2-dependent E3 ubiquitin-protein ligase. May be involved in endocytic trafficking.

## Background

Acts as an E2-dependent E3 ubiquitin-protein ligase. May be involved in endocytic trafficking.

---

Note: This product is for in vitro research use only