

TMEM43 Rabbit mAb

Catalog No: #60211

Package Size: #60211-1 50ul #60211-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TMEM43 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	TMEM43 Antibody detects endogenous levels of total TMEM43
Immunogen Description	A synthesized peptide derived from human TMEM43
Other Names	ARVC5; ARVD5; EDMD7; LUMA; Tmem43; Transmembrane protein 43;
Accession No.	Uniprot:Q9BTV4
Uniprot	Q9BTV4
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of TMEM43 expression in Jurkat cell lysate.

Product Description

May have an important role in maintaining nuclear envelope structure by organizing protein complexes at the inner nuclear membrane. Required for retaining emerin at the inner nuclear membrane.

Background

May have an important role in maintaining nuclear envelope structure by organizing protein complexes at the inner nuclear membrane. Required for retaining emerin at the inner nuclear membrane.

Note: This product is for in vitro research use only