

UCHL5IP Rabbit mAb

Catalog No: #60219

Package Size: #60219-1 50ul #60219-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	UCHL5IP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Human
Specificity	UCHL5IP Antibody detects endogenous levels of total UCHL5IP
Immunogen Description	A synthesized peptide derived from human UCHL5IP
Other Names	Haus7; HSXQ28ORF; UCHL5IP; UIP1;
Accession No.	Uniprot:Q99871
Uniprot	Q99871
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 IP 1:50

Images



Western blot analysis of UCHL5IP expression in Jurkat cell lysate.

Product Description

Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.

Background

Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.

Note: This product is for in vitro research use only