

Annexin A7 Rabbit mAb

Catalog No: #60223

Package Size: #60223-1 50ul #60223-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Annexin A7 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	Annexin A7 Antibody detects endogenous levels of total Annexin A7
Immunogen Description	A synthesized peptide derived from human Annexin A7
Other Names	Annexin A7; Annexin VII; ANX7; ANXA7; SNX; Synexin;
Accession No.	Uniprot:P20073
Uniprot	P20073
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of Annexin A7 expression in (1) Jurkat cell lysate; (2) Raw264.7 cell lysate.

Product Description

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Background

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Note: This product is for in vitro research use only