

## PI3 Kinase p110 beta Rabbit mAb

Catalog No: #60224

Package Size: #60224-1 50ul #60224-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PI3 Kinase p110 beta Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	PI3 Kinase p110 beta Antibody detects endogenous levels of total PI3 Kinase p110 beta
Immunogen Description	A synthesized peptide derived from human PI3 Kinase p110 beta
Other Names	p110Beta; PI3K; PI3K beta; PI3KCB; PIK3C1; Pik3cb;
Accession No.	Uniprot:P42338
Uniprot	P42338
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of PI3 Kinase p110 beta expression in (1) 293 cell lysate; (2) Mouse Brain lysate.

## Product Description

Phosphorylates PtdIns, PtdIns4P and PtdIns(4,5)P2 with a preference for PtdIns(4,5)P2.

## Background

Phosphorylates PtdIns, PtdIns4P and PtdIns(4,5)P2 with a preference for PtdIns(4,5)P2.

---

Note: This product is for in vitro research use only