

## SNX4 Rabbit mAb

Catalog No: #60225

Package Size: #60225-1 50ul #60225-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	SNX4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	SNX4 Antibody detects endogenous levels of total SNX4
Immunogen Description	A synthesized peptide derived from human SNX4
Other Names	ATG24B; snx4; Sorting nexin 4;
Accession No.	Uniprot:O95219
Uniprot	O95219
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000 IP 1:50

## Images



Western blot analysis of SNX4 expression in (1) A431 cell lysate; (2) Raw264.7 cell lysate.

## Product Description

May be involved in several stages of intracellular trafficking. Plays a role in recycling endocytosed transferrin receptor and prevent its degradation.

## Background

May be involved in several stages of intracellular trafficking. Plays a role in recycling endocytosed transferrin receptor and prevent its degradation.

---

Note: This product is for in vitro research use only