

ACTR1B Rabbit mAb

Catalog No: #60230

Package Size: #60230-1 50ul #60230-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ACTR1B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human
Specificity	ACTR1B Antibody detects endogenous levels of total ACTR1B
Immunogen Description	A synthesized peptide derived from human ACTR1B
Other Names	ACTR1B; ARP1B; Beta centractin; CTRN2; PC3;
Accession No.	Uniprot:P42025
Uniprot	P42025
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of ACTR1B expression in 293 cell lysate.

Product Description

Component of a multi-subunit complex involved in microtubule based vesicle motility. It is associated with the centrosome.

Background

Component of a multi-subunit complex involved in microtubule based vesicle motility. It is associated with the centrosome.

Note: This product is for in vitro research use only