

PRPF4 Rabbit mAb

Catalog No: #60231

Package Size: #60231-1 50ul #60231-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PRPF4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human Mouse Rat
Specificity	PRPF4 Antibody detects endogenous levels of total PRPF4
Immunogen Description	A synthesized peptide derived from human PRPF4
Other Names	hPrp4; HPRP4P; PRP4; Prp4p; PRPF4; RP70; SNRNP60;
Accession No.	Uniprot:O43172
Uniprot	O43172
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Images



Western blot analysis of PRPF4 expression in (1) Raji cell lysate; (2) RAW 264.7 cell lysate; (3) C6 cell lysate.

Product Description

Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.

Background

Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.

Note: This product is for in vitro research use only