

PMF1 Rabbit mAb

Catalog No: #60233

Package Size: #60233-1 50ul #60233-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PMF1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human Mouse Rat
Specificity	PMF1 Antibody detects endogenous levels of total PMF1
Immunogen Description	A synthesized peptide derived from human PMF1
Other Names	Pmf1; Polyamine modulated factor 1;
Accession No.	Uniprot:Q6P1K2
Uniprot	Q6P1K2
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IP 1:50

Images



Western blot analysis of PMF1 expression in (1) HeLa cell lysate; (2) RAW 264.7 cell lysate.

Product Description

Part of the MIS12 complex which is required for normal chromosome alignment and segregation and kinetochore formation during mitosis. May act as a cotranscription partner of NFE2L2 involved in regulation of polyamine-induced transcription of SSAT.

Background

Part of the MIS12 complex which is required for normal chromosome alignment and segregation and kinetochores formation during mitosis. May act as a cotranscription partner of NFE2L2 involved in regulation of polyamine-induced transcription of SSAT.

Note: This product is for in vitro research use only