

RPS8 Rabbit mAb

Catalog No: #60242

Package Size: #60242-1 50ul #60242-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RPS8 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	RPS8 Antibody detects endogenous levels of total RPS8
Immunogen Description	A synthesized peptide derived from human RPS8
Other Names	RPS8;
Accession No.	Uniprot:P62241
Uniprot	P62241
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 IP 1:50 FC 1:50

Images



Western blot analysis of RPS8 expression in (1) 293 cell lysate; (2) RAW 264.7 cell lysate.

Product Description

Belongs to the ribosomal protein S8e family.

Background

Belongs to the ribosomal protein S8e family.

Note: This product is for in vitro research use only