

TCTEX1 Rabbit mAb

Catalog No: #60254

Package Size: #60254-1 50ul #60254-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TCTEX1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse
Specificity	TCTEX1 Antibody detects endogenous levels of total TCTEX1
Immunogen Description	A synthesized peptide derived from human TCTEX1
Other Names	AGS2; DYNLT1; TCTEL1; Tctex1;
Accession No.	Uniprot:P63172
Uniprot	P63172
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Images



Western blot analysis of TCTEX1 expression in HeLa cell lysate.

Product Description

Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function.

Background

Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function.

Note: This product is for in vitro research use only