

## TCEB2 Rabbit mAb

Catalog No: #60259

Package Size: #60259-1 50ul #60259-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TCEB2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	TCEB2 Antibody detects endogenous levels of total TCEB2
Immunogen Description	A synthesized peptide derived from human TCEB2
Other Names	EloB; Elongin B; SIII; SIII p18; Tceb2;
Accession No.	Uniprot:Q15370
Uniprot	Q15370
Calculated MW	16kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details


WB 1:500~1:2000

## Images



Western blot analysis of TCEB2 expression in Jurkat cell lysate.

All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



## Product Description

SIII, also known as elongin, is a general transcription elongation factor that increases the RNA polymerase II transcription elongation past template-encoded arresting sites.

## Background

SIII, also known as elongin, is a general transcription elongation factor that increases the RNA polymerase II transcription elongation past template-encoded arresting sites.

Note: This product is for in vitro research use only