

ABCD1 / ALD Rabbit mAb

Catalog No: #60268

Package Size: #60268-1 50ul #60268-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ABCD1 / ALD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human Mouse Rat
Specificity	ABCD1 / ALD Antibody detects endogenous levels of total ABCD1 / ALD
Immunogen Description	A synthesized peptide derived from human ABCD1 / ALD
Other Names	ABC42; Abcd1; ALD; Aldgh; ALDP; AMN;
Accession No.	Uniprot:P33897
Uniprot	P33897
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200 FC 1:50

Images



Western blot analysis of ABCD1 / ALD in HepG2 cell lysate.

Product Description

Probable transporter. The nucleotide-binding fold acts as an ATP-binding subunit with ATPase activity.

Background

Probable transporter. The nucleotide-binding fold acts as an ATP-binding subunit with ATPase activity.

Note: This product is for in vitro research use only