

## Cathepsin H Rabbit mAb

Catalog No: #60273

Package Size: #60273-1 50ul #60273-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |  |
|-----------------------|--|
| Product Name          | Cathepsin H Rabbit mAb   |
| Host Species          | Rabbit   |
| Clonality             | Monoclonal   |
| Isotype               | Rabbit IgG   |
| Purification          | Affinity-chromatography  |
| Applications          | WB IHC ICC/IF  |
| Species Reactivity    | Human  |
| Specificity           | Cathepsin H Antibody detects endogenous levels of total Cathepsin H                                |
| Immunogen Description | A synthesized peptide derived from human Cathepsin H   |
| Other Names           | ACC4; ACC5; cathepsin B3; cathepsin BA; CPSB; CTSH;  |
| Accession No.         | Uniprot:P09668   |
| Uniprot               | P09668   |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.                     |

## Application Details

WB 1:500~1:2000 IHC 1:50~1:200 ICC/IF 1:50~1:200

## Images



Western blot analysis of Cathepsin H in HepG2 cell lysate.

## Product Description

Important for the overall degradation of proteins in lysosomes.

## Background

Important for the overall degradation of proteins in lysosomes.

---

Note: This product is for in vitro research use only