

## IL22 Rabbit mAb

Catalog No: #60276

Package Size: #60276-1 50ul #60276-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	IL22 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	IL22 Antibody detects endogenous levels of total IL22
Immunogen Description	A synthesized peptide derived from human IL22
Other Names	IL21; IL22; ILD110; ILTIF; interleukin 21; Interleukin 22; TIFa; TIFIL23; zcyto18;
Accession No.	Uniprot:Q9GZX6
Uniprot	Q9GZX6
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of IL22 in Recombinant Human IL22 protein cell lysate.

## Product Description

Cytokine that contributes to the inflammatory response in vivo.

## Background

Cytokine that contributes to the inflammatory response in vivo.

---

Note: This product is for in vitro research use only