

IL9 Rabbit mAb

Catalog No: #60277

Package Size: #60277-1 50ul #60277-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IL9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	IL9 Antibody detects endogenous levels of total IL9
Immunogen Description	A synthesized peptide derived from human IL9
Other Names	HP40; IL9; Interleukin 9 MCGF; P40; TCGF 3;
Accession No.	Uniprot:P15248
Uniprot	P15248
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of IL9 in Recombinant Human IL9 protein cell lysate.

Product Description

Supports IL-2 independent and IL-4 independent growth of helper T-cells.

Background

Supports IL-2 independent and IL-4 independent growth of helper T-cells.

Note: This product is for in vitro research use only