

## IL26 Rabbit mAb

Catalog No: #60278

Package Size: #60278-1 50ul #60278-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	IL26 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	IL26 Antibody detects endogenous levels of total IL26
Immunogen Description	A synthesized peptide derived from human IL26
Other Names	IL26; Interleukin 26;
Accession No.	Uniprot:Q9NPH9
Uniprot	Q9NPH9
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Application Details

WB 1:500~1:2000

## Images



Western blot analysis of IL26 in His-tagged IL26 cell lysate.

## Product Description

May play a role in local mechanisms of mucosal immunity and seems to have a proinflammatory function. May play a role in inflammatory bowel disease.

## Background

---

May play a role in local mechanisms of mucosal immunity and seems to have a proinflammatory function. May play a role in inflammatory bowel disease.

---

Note: This product is for in vitro research use only