

CD125 Rabbit mAb

Catalog No: #60281

Package Size: #60281-1 50ul #60281-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD125 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	CD125 Antibody detects endogenous levels of total CD125
Immunogen Description	A synthesized peptide derived from human CD125
Other Names	CD125; CDw125; HSIL5R3; IL5R; IL5RA;
Accession No.	Uniprot:Q01344
Uniprot	Q01344
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000

Images



Western blot analysis of CD125 in Human serum cell lysate.

Product Description

This is the receptor for interleukin-5. The alpha chain binds to IL5.

Background

This is the receptor for interleukin-5. The alpha chain binds to IL5.

Note: This product is for in vitro research use only