

Dopamine beta Hydroxylase Rabbit mAb

Catalog No: #60283

Package Size: #60283-1 50ul #60283-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Dopamine beta Hydroxylase Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	Dopamine beta Hydroxylase Antibody detects endogenous levels of total Dopamine beta Hydroxylase
Immunogen Description	A synthesized peptide derived from human Dopamine beta Hydroxylase
Other Names	dbh; DBM; Dopamine beta hydroxylase; Dopamine beta monooxygenase; Soluble dopamine beta-hydroxylase;
Accession No.	Uniprot:P09172
Uniprot	P09172
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Application Details

WB 1:500~1:2000 IHC 1:50~1:100

Images



Western blot analysis of Dopamine beta Hydroxylase expression in SH-SY5Y cell lysate.

Product Description

Conversion of dopamine to noradrenaline.

Background

Conversion of dopamine to noradrenaline.

Note: This product is for in vitro research use only